



Focus on
**Clinical
Consumables**



SciLabware

AZLON® MBL® PYREX® QUICKFIT® WHEATON®

Clinical Consumables

Introduction

SciLabware is one of the world's leading manufacturer's of laboratory glass and plasticware, supplying the scientific community with world renown brands including Azlon®, Pyrex®, Quickfit®, Wheaton® and Ramboldi®.

SciLabware offers a range of products specifically designed for the clinical laboratory, offering a cost effective choice of sample containers, blood collection tubes, inoculation needles and loops, transfer pipettes, pipette tips, tubes, gallipots, wash bottles, wash bottles and media vials and bottles to hospital pathology laboratories.



Gallipots

Azlon®

- Medical grade polypropylene.
- Suitable for sampling, short term storage, micro-analysis and many other applications.
- Moulded graduations to capacity.
- Available sterile and non-sterile.



Wash Bottles

Azlon®

- Pre-labelled with gram stain solutions (Crystal violet, Lugol's, Gram iodine, Gram colour/de-colourising and Safranin) and bleach.
- Patented Driplok® vapour venting valve prevents solvent drips.
- Includes new GHS symbols, NFPA diamond, risk phrases and safety information.



Sample Collection Containers

RAMBOLDI

- Manufactured from virgin, medical grade, polystyrene or polypropylene.
- Ideal for sample containment and storage.
- Supplied sterile for in vitro use only.
- CE marked in accordance with European Directive 98/79/EC
- All containers are leak tested.
- Available in 7ml, 30ml 60ml, 100ml 150ml and 250ml sizes.
- 7ml and 30ml available with or without boric acid.
- All containers available with no, plain, or printed label.
- Packed in inner bags for easy distribution.

Primary Urine Tubes



- Manufactured from virgin, medical grade, polystyrene.
- Ideal for use with automated urine analysis machines.
- Removes the need for unnecessary decanting of sample in the laboratory.
- Supplied sterile for in vitro use only.
- CE marked in accordance with European Directive 98/79/EC.
- All tubes are leak tested.
- Option of plain tube, or pre-dosed with boric acid.
- Printed label for patients details.
- Packed either in inner bags of 50 for easy distribution, or individually wrapped.
- Collection cup supplied with all tube options.



Plastic Pasteur Pipettes



- Manufactured from non-toxic LDPE.
- Ideal for safe and easy transfer of liquids and samples.
- Available in 1ml, 1ml fine tip, 3ml and 10ml sizes.
- Graduations on 1ml and 3ml for easy measurement.
- Supplied non-sterile or sterile.
- Sterile pipettes available in inner packs, individually wrapped, in sets of 10 or 20.



Autoclave and Specimen Bags



Printed Autoclave Bags

- High Temperature autoclave bags, suitable for use up to 135°C.
- Manufactured from heavy gauge polypropylene.
- Available in 3 sizes.

Specimen Bags

- Grip seal polyethylene bags suitable for sampling purposes.
- Available plain or with biohazard sign for hazardous samples.
- All specimen bags have document wallets as standard.





SciLabware

AZLON® MBL® PYREX® QUICKFIT® WHEATON®

For more information visit our website
www.scilabware.com

Email us with your questions
info@scilabware.com

SciLabware Limited Unit 4, Riverside 2,
Campbell Road, Stoke on Trent, ST4 4RJ, UK.
Tel: +44 (0)1782 444406 Fax: +44 (0)1782 940436



A DWK LIFE SCIENCES COMPANY