

Centrifuges made by Funke-Gerber
approved – exact – reliable

**FUNKE
GERBER**

SuperVario N

Multi-Purpose Centrifuge for the Dairy Industry



Multi-Purpose Centrifuge SuperVario N

Gerber-Centrifuge up to 36 Butyrometers, Babcock fat determination, Röse-Gottlieb fat determination* and solubility index determination (ATPI)

**Pay attention to existing safety rules!*

Features:

- **Easy Handling:**
 - 5 selectable modes for centrifugation** (Gerber test, Röse-Gottlieb, Babcock, Solubility and Free Mode)
 - Free Set-up (Free Mode) of**
 1. **Speed/RCA**** (adjustable from 600 to 1200 rpm, the corresponding RCA is shown simultaneously)
 - **RCA = Relative Centrifugal Acceleration
 2. **Temperature** (adjustable from amb. temp. to 68°C)
 3. **Centrifugation time** (adjustable from 1 to 99 minutes)
- **Centrifugation drum of stainless steel**
- **Security lid locking device**
- **Automatic unbalance switch off**
- **Electrical brake**

Technical Specification:



Connection: 230V/50...60 Hz/1200 VA
 Weight (empty): 26 kg
 Full height: 460 mm
 Operation height: 370 mm
 Speed range: 600...1200 UpM
 Temperature range: Amb. temperature to 68°C



Order Data:

Name	Art. No.
Centrifuge SuperVario-N (without heads)	3680

Accessories:

Name		Art. No.
A 	Head A (up to 36 Butyrometer/up to 18 Babcock–bottles)	3685
	Butyrometer–bucket of cast aluminium	3631
	Babcock–bucket	3632
B 	Head B (Protective drum for 8 Mojonnier tubes)	3686
C 	Head C (up to 6 solubility index tubes)	3687
	Solubility index tube bucket	3633