







Scharlau filter tips

RNase & DNase free · Pyrogen free · Autoclavable at 121°C/15min
Graduation marks · High grade of transparency

Scharlau filter tips compatibility table*



					
00PF0110-2 0,1-10µl Rack 10x96	PFXL0110-2 0,1-10µl Rack 10x96	000PF120-2 1-20µl Rack 10x96	00PF1100-2 1-100µl Rack 10x96	00PF1200-2 1-200µl Rack 10x96	PF001000-2 100-1000µl Rack 8x96
00PF0110-1 0,1-10µl Bag 1.000 units	PFXL0110-1 0,1-10µl Bag 1.000 units	000PF120-1 1-20µl Bag 1.000 units	00PF1100-1 1-100µl Bag 1.000 units	00PF1200-1 1-200µl Bag 1.000 units	PF10010001 100-1.000µl Bag 1.000 units

Model	Trademark	10µl	10µl XL	20µl	100µl	200µl	1000µl
Transferpette®	Brand	✓	✓			✓	✓
CAPP Bravo	CAPP	✓	✓	✓	✓	✓	✓
CAPP Aero		✓	✓	✓	✓	✓	✓
Research® Plus	Eppendorf	✓	✓	✓	✓	✓	✓
Finnpipette Digital	Finnpipette	✓	✓		✓	✓	✓
Pipetman Classic™	Gilson	✓	✓	✓	✓	✓	✓
Rainin Classic	Rainin	✓	✓		✓	✓	✓
Proline® Plus	Sartorius	✓	✓		✓	✓	✓
mLINE		✓	✓		✓	✓	✓
Automatic pipettes	VWR	✓	✓		✓	✓	✓
	Labnet	✓	✓		✓	✓	✓
	HTL	✓	✓		✓	✓	✓

*Compatibilities are mere suggestions of the supplier.