

# Scharlau heating blocks



- Made of an anodized aluminium alloy
- Suitable for heating round bottom flasks from 10ml to 1L on a stirring hotplate
- Compatible with any stirring hotplate and can be placed anywhere on top of it
- Lateral hole allows block temperature measurement using a PT100 probe
- Block bottom holes allow the release of fumes in case liquid should accidentally leak from the flask
- Replace oil baths and heating mantles increasing speed, cleanness and safety



*The wise choice*



Hotplate not included

## Ordering information

ART. NO.	DESCRIPTION
CAL-000010	Aluminium heating block for 10ml flask
CAL-000025	Aluminium heating block for 25ml flask
CAL-000050	Aluminium heating block for 50ml flask
CAL-000100	Aluminium heating block for 100ml flask
CAL-000250	Aluminium heating block for 250ml flask
CAL-000500	Aluminium heating block for 500ml flask
CAL-001000	Aluminium heating block for 1000ml flask

