

# PRODUCT INFO

Scharlau  
by  
AFORA

Scharlab  
THE LAB SOURCING GROUP

CON  
chemicals ltd

## Accreditation for volumetric glassware according to ISO 17025:2017

Our glassware workshop, Afora, is equipped with its own Calibration Laboratory in its facilities, accredited by ENAC\* in accordance with the criteria included in the UNE-EN ISO/IEC 17025 standard in its current version.



certiglass



> **SCIENTIFIC GLASS FACTORY WITH ITS OWN ACCREDITED VOLUMETRIC CALIBRATION LABORATORY ACCORDING TO ISO 17025:2017**

> **SERVICE FOR CALIBRATION LABORATORIES:**

- IMPROVED CONTROL
- IMPROVED RECOGNITION
- IMPROVED EFFICIENCY

> **DISCOVER THE MANUFACTURING AND CERTIFICATION PROCESS IN OUR INSTALATIONS:**



Afora Volumetric Calibration Laboratory is fully equipped with all necessary means, to both guarantees controlled environmental conditions as well as accurate and reliable measurements, which metrologically grants all the results reported in our certificates. We have two precision automatic machines manufactured and tuned by the German Fraunhofer Institute, for the calibration of volumetric flasks (petrol station flasks according to Spanish regulations), pipettes and burettes.

Regardless of the Scope of Accreditation for which we are audited by ENAC, we can issue certificates for batch volumetric pieces, individual certificates and certificates according to the requirements of USP (The United States Pharmacopeia).



PI-ENENG22



\*ENAC is a signatory of all international agreements of the EA, ILAC and IAF. This means that a report or certificate issued under ENAC accreditation will be recognised by the other signatories across the world.